UNMH: Validation Tool: Elastomeric Masks

nstructions: This validation tool must be completed in order for ind Validation Key: $\mathbf{O} = \text{Observed practice}$. $\mathbf{S} = \text{Simulated demonstrates}$			= Verbalizes und	
Criteria	Compo	etent No	Verification Method	
 Locate clean masks and filters in SPD containers. Check out mask. 			V	
 Inspect mask Check to ensure filter is not soiled or damaged – if so, discard and obtain new filter Attach both filters with a quarter turn 			O, S	
 Appropriately don mask Appropriately conduct positive and negative pressure seal check 			O, S	
 Secure a surgical mask over the Respirator exhalation port Appropriately Doff mask Place respirator interior side down in designated bucket 			O, S S	
 Don new gloves and surgical mask if transporting patient (to PACU etc) 			S	
 Remove and dispose of surgical mask from EMR while wearing gloves. Perform hygiene of gloved hands. Clean exterior of mask with 1st wipe being careful not to wipe inside cartridge or saturate the filter material. 			S	
 If elastomeric has been doffed for the last time in the dayremove filters and continue wiping all outside surfaces of Filter. (If not, then keep filters on.) Place on clean chux 			S	
 With 2nd wipe, clean inside of mask and wipe down straps Place in level 2 container and take to SPD for additional cleaning/disinfection. 			S	
 Document mask checked back in Return filters to storage area have completed these requirements and understand the information p 			V	
have completed these requirements and understand the information produced in the infor	oresented r	egardii	ng Elastomeric do	nning, doffin
Employee signature		Dat	e	
Batcave/Unit Based Educator/Unit Director/Supervisor				
Printed Name:		Date		

Plan: This section to be completed <u>only</u> if an improvement plan is needed. Needs additional supervised experience Y/N Limit scope of practice

Needs additional supervised experience Y/N Needs further review Y/N

Limit scope of practice Y/N Initiate disciplinary action Y/N